

NRAS, WT1, NPM1 and IDH mutation detection for molecular MRD using Flow Cytometry in pediatric AML

Margarita Maurer-Granofszky^{1,2}, Stefan Köhrer², Astrid Mecklenbraeuer², Linus Bosaeus³, Lei Chen³, Michael Dworzak^{1,2,4}

1. St. Anna Children's Cancer Research Institute (CCRI), Vienna, Austria
2. Labdia Labordiagnostik, Vienna, Austria
3. Rarity Bioscience, Uppsala, Sweden
4. St. Anna Children's Hospital, Department of Pediatrics and Adolescent Medicine, Medical University of Vienna, Vienna, Austria




INTRODUCTION

MRD is a key prognostic indicator in pediatric leukemia. Molecular methods complement multiparametric Flow Cytometry (MFC), enhancing sensitivity and early relapse detection. We present data from a novel molecular assay using Flow Cytometry for the multiplex detection of driver mutations, including NRAS, WT1, NPM1, and IDH1/2.

METHOD

In this retrospective non-interventional study, the initial diagnostic work-up contained mutation screening using targeted Next Generation Sequencing (Trusight Myeloid Panel, Illumina, CA, USA) at the BFM reference laboratory for molecular diagnostics in pediatric acute myeloid leukemia in Essen, Germany. Follow-up samples had previously been analyzed using MFC (DxFLEX, Beckman Coulter, CA, USA) and, when applicable, DNA-based genomic breakpoint specific qPCR following EuroMRD guidelines. For this study, patient samples from diagnosis and follow-up (including relapse timepoints) preserved in liquid nitrogen were analyzed and compared using a new molecular MRD method superRCA (Rarity Bioscience, Uppsala, Sweden) for potential driver gene mutations identified at diagnosis. Analysis was done over two days by two operators, one previously untrained.



RESULTS

- 23 bone marrow samples from 5 pediatric AML patients analyzed for potential driver mutations in NRAS (4), WT1 (1), NPM1 (2), IDH1 (1)
- superRCA assays mean Wild-type control= 2,11E-05 (0-8,16E-05)
- superRCA between operators CV= 4.11% (0.33-8.68%)
- Correlation between superRCA and NGS, mean variation of 3,55% across all targets, Pearson correlation of $r=0.93$, $p=0,0009$




Figure 1. Mean Wild type background, control samples in replicates




Figure 3: Correlation between MFC and superRCA at diagnosis, relapse, and during follow-up in a patient (Vie3) who developed treatment resistance after relapse. Divergence between MFC and superRCA suggests the presence of subclonal mutations. Abbreviations: IND1, post-Induction 1; IND2, post-Induction 2; CON1, post-Consolidation 1; Rel, Relapse

Table 1. Sample analysis including MFC, PCR, and mutation analysis superRCA (NGS)

Patient		Diag.	IND1	IND2	CON1	Relapse	Rel_IND1	Rel_IND2
Vie 1	MFC	38%	0	-	-	-	-	-
	NPM1 Mut A	28,28% (28,84%)	0,06%	-	-	-	-	-
	WT1	21,00% (24,22%)	0,03%	-	-	-	-	-
	c.1078_1091dup							
Vie 2	MFC	34%	0%	-	-	-	-	-
	NPM1 Mut B	47,54% (37,7%)	0,10%	-	-	-	-	-
Vie 3	MFC	61%	-	-	-	6,30%	20,98%	33,60%
	NRAS p.G12S	35,88% (38,2%)	-	-	-	0,207%	2,63%	0,143%
Vie 4	MFC	95%	-	0%	0%	-	-	-
	IDH1 p.R132C	43,96% (45,5%)	-	0,09%	0,39%	-	-	-
	NRAS p.G12D	30,39% (34,8%)	-	0,02%	0,08%	-	-	-
	NRAS p.Q61K	8,35% (8,2%)	-	0,00%	0,01%	-	-	-
Vie 5	MFC	84%	0%	0%	0%	-	-	-
	PCR - MLL-MLLT3	-	Neg	Neg	Neg	-	-	-
	NRAS p.G12A	40,5% (46,85%)	0,00%	0,00%	0,00%	-	-	-




Figure 2. Correlation between NGS and superRCA MAF all targets




Figure 4. Correlation between MFC and superRCA at diagnosis and two follow-up time points of Patient Vie4. Instances where MFC is negative but superRCA remains positive may indicate either that blast counts are below MFC's detection threshold or that mutations exist in mature cells, which are not captured when monitoring MRD via MFC. Abbreviations: IND1, post-Induction 1; IND2, post-Induction 2; CON1, post-Consolidation 1.

CONCLUSION

Analysis of molecular MRD using flow cytometry is viable approach and may particularly add value to pediatric leukemia patients with identified driver mutations. The superRCA assay shows strong correlation to the NGS data at diagnosis and demonstrated a very low limit of detection for molecular MRD analysis at follow up. This warrants further evaluation in prospective studies.